


## Putting your earmould in



1. Please stand in front of the mirror as it may be easier to put the earmould into your ear.


2. Hold the earmould by the handle using your thumb and forefinger (the part where it looks like a tea cup handle).

Tea Cup Handle



3. Take the earmould up to your ear and place the top end of the earmould into the fold of your ear.

Top End



4. Place the bottom end of the earmould into your ear canal .



**Bottom end**

5. If you are having difficulty inserting the earmould into your ear canal, use your other hand to gently pull the earlobe downwards as this opens the ear canal.



6. Press the earmould firmly into your ear.



7. Place the hearing aid behind you ear. The tube and hearing aid should be sitting in your ear as demonstrated in the picture. If possible ask someone to check if the earmould and hearing aid is sitting in your ear as shown in the picture.

